

**MWR-fix External fixator**

### **MWR-Fix**

**MWR-fix is an External fixator for installed bone fractures and bone lengthening.**

**Complete. 10.000.10**

**It's an External fixator device -with 2 fixing units- allows to the surgeon the possibility of installing bone fractures and using the fixing bars at the level of two dimensions up to 360° of each.**


**Content of**

- **Fixing Unite, 2 pcs , The possibility of the use of fixing bars from 1,5 mm to 7mm diameter with the compress. 10.100.10**



**-Fixing Joins , 2pcs, The possibility of the use of fixing bars at the level of two dimensions up to 360° of each. 10.200.10**





**MWR-fix External fixator**

**-Fixing Transmisson, 2PCS,  
10.300.10**




**Compressor and -Lengthening Axis, Lengthening from 1 – 12 CM  
10.400.10**



**Ruler Axis,  
10.500.10**





MWRS-fix External fixator

MWRS-fix is an External fixator for installed bone fractures and bone lengthening.  
Complete. 10.010.10

It's an External fixator device -with 4 fixing units- allows to the surgeon the possibility of installing bone fractures and using the fixing bars at the level of three dimensions up to 360° of each.


Content of

- Fixing Unite, 4 pcs , The possibility of the use of fixing bars from 1,5 mm to 7mm diameter with the compress . 10.110.10



-Fixing Joins , 4pcs, The possibility of the use of fixing bars at the level of two dimensions up to 360° of each. 10.210.10





**MWRS-fix External fixator**

**Joins Holder,2PCS.**


**10.310.10**



**-Fixing Transmisson, 2PCS,**

**10.410.10**






**MWRS-fix External fixator**

**Compressor and -Lengthening Axis, Lengthening from 1 – 12 CM  
10.510.10**






**Ruler Axis,  
10.610.10**






MWRS-fix External fixator






MWRS-fix External fixator






MWRS-fix External fixator







MWRS-FIX MINI External fixator

### MWR-FIX MINI

MWR-fix Mini is an External fixator for installed bone fractures and bone lengthening.

Complete. 10.020.10

It's an External fixator device -with 4 fixing units- allows to the surgeon the possibility of installing bone fractures and using the fixing bars at the level of three dimensions up to 360° of each.


Content of

- Fixing Unite, 4 pcs , The possibility of the use of fixing bars from 1 mm to 4 mm diameter. with the compress 10.120.10



-Fixing Joins , 4pcs, The possibility of the use of fixing bars at the level of two dimensions up to 360° of each. 10.220.10





**MWRS-FIX MINI External fixator**

**Joins Holder,2PCS.**


**10.320.10**



**-Fixing Transmisson, 2PCS,**

**10.420.10**






**MWRS-FIX MINI External fixator**

**Compressor and -Lengthening Axis, Lengthening from 1 – 12 CM  
10.520.10**







**Ruler Axis,  
10.620.10**







MWRS-FIX MINI External fixator





MWRS-FIX MINI External fixator





GH-MED  
Hannaström Medicinska instrument HB  
Adress: Holbergsgatan.51 , 16849 Bromma ,Stockholm Sweden  
Tele/Fax:+46(0)8-579 72 900 Webside.[www.gh-med.se](http://www.gh-med.se) E-mail:[info@gh-med.se](mailto:info@gh-med.se)